


# ShowCase®

Clinical Case Presentation Tools

DICOM Viewer Products

Specifications & Prices  
**SHOWCASE 5**

## ShowCase Premier



### Features

- Display DICOM images from DICOM CD/DVD/disk directories or individual DICOM files.
- **Basic image measurement tool set – line, ellipse and free form areas; Doppler and M-Mode line tools, Spectral Doppler trace tool.**
- Display measurements /clinical data transferred from ultrasound machines as DICOM Structured Reports.
- **Export studies to CD or DVD with Onboard Viewer**
- Margin clips or crop images for export to slides.
- User selectable sizes and compare modes.
- Multi-frame image "play" and "step" modes with speed control for "play"
- Study de-identification tools for research, presentations and teaching files, accreditation examples.
- Direct export to PowerPoint® presentations
- Contrast / gamma adjust images
- Export images in bmp, tiff, jpeg, avi formats
- Copy/paste/delete images to make new custom series for teaching, marketing, special review.
- Save As and Copy to DICOM media
- Print reference images to desktop printers
- Annotations and erase image areas
- **CARDIOLOGY** - Stress Echo displays, including export stage or view layouts to AVI files. IVUS displays

### Option -- ShowCase Connect (DICOM Networking)



The Connect module lets you to store images and DICOM structured reports transferred over a network from ultrasound machines (a DICOM image store SCP), then view them with Premier

You can also Query/Retrieve studies from qualified DICOM servers (PACS servers)

### Option-- Portable USB Hardware License Key



A hardware license key provides easy licensing without an Internet connection. It is also handy for installing ShowCase on multiple machines and then using one at a time – a portable license.

Product	Part #	Price
ShowCase Premier	PN050-e	900
ShowCase Premier – USB Portable License	PN050-usb	960
ShowCase Premier with Connect	PC050-e	1500

## ShowCase Basic



### Features

- Display DICOM images from DICOM CD/DVD/disk directories or individual DICOM files.
- Display measurements /clinical data transferred from ultrasound machines as DICOM Structured Reports.
- Margin clips or crop images for export to slides.
- User selectable sizes and compare modes.
- Multi-frame image "play" and "step" modes with speed control for "play"
- Study de-identification tools for research, presentations and teaching files, accreditation examples.
- Direct export to PowerPoint® presentations
- Contrast / gamma adjust images
- Export images in bmp, tiff, jpeg, avi formats
- Copy/paste/delete images to make new custom series for teaching, marketing, special review.
- Save As and Copy to DICOM media
- Print reference images to desktop printers
- Annotations and erase image areas
- **CARDIOLOGY** - Stress Echo displays, including export stage or view layouts to AVI files. IVUS displays

### Option -- ShowCase Connect (DICOM Networking)



The Connect module allows you to store images and DICOM structured reports transferred over a network from ultrasound machines (a DICOM image store SCP), then view them with ShowCase

You can also Query/Retrieve studies from qualified DICOM servers (PACS servers)



### Option-- Portable USB Hardware License Key



A hardware license key provides easy licensing without an Internet connection. It is also handy for installing ShowCase on multiple machines and then using one at a time – a portable license.

Product	Part #	Price
ShowCase Basic	SN050-e	500
ShowCase Basic – USB Portable License	SN050-usb	560
ShowCase Basic with Connect	SC050-e	1100

## Additional Products

	Product	Part #	Price
	<b>Backup CD and a hardcopy User Guide –</b> (for any product). Note that USB licenses include installation software on the USB drive. Non-USB licenses (Internet or manually activated) are priced for delivery via electronic software download.	GCD050	20
	<b>Additional User Guides</b>	SCGD05	12